

MXD4-MAX dimerization protein 4 Gene

Reference: PTXBC002713

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAD4; MST149; MSTP149; bHLHc12; max dimerization protein 4; Mad4 homolog; class C basic helix-loop-helix protein 12; max dimerizer 4; max-associated protein 4; max-interacting transcriptional repressor MAD4		

For research use only.

GENE DESCRIPTION	MAX dimerization protein 4
GENE ID	10608
ORIGIN SPECIES	Human
ACCESSIONS	BC002713
VECTOR	pDONR223
NCBI SYMBOL	MXD4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctgaactccctgctgatctgctgagggcgccgagtagctggagcgcagggatcgagagccgagcagcgctacgctcggtgctgccctcgacggcgactcgccagggagagaaacaaaggcgccgctggtgcgaaggccccgaacaacaggcctctcagtga