



Reference: PTXBC002730

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DD83.1; RDH17; Rsdr1; SDR1; SDR16C1; retSDR1; short-chain dehydrogenase/reductase 3; dehydrogenase/reductase (SDR family) member 3; dehydrogenase/reductase member 3; retinal short-chain dehydrogenase/reductase 1; retinol dehydrogenase 17; short chain dehy		

For research use only.

GENE DESCRIPTION	dehydrogenase/reductase (SDR family) member 3
GENE ID	9249
ORIGIN SPECIES	Human
ACCESSIONS	BC002730
VECTOR	pDONR223
NCBI SYMBOL	DHRS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1