

CRIP1-cysteine-rich protein 1 (intestinal) Gene

Reference: PTXBC002738

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRHP; CRIP; CRP-1; CRP1; cysteine-rich protein 1; cysteine-rich heart protein; cysteine-rich intestinal protein; cysteine-rich protein 1 (intestinal); cysteine rich protein 1		

For research use only.

GENE DESCRIPTION	cysteine-rich protein 1 (intestinal)
GENE ID	1396
ORIGIN SPECIES	Human
ACCESSIONS	BC002738
VECTOR	pDONR223
NCBI SYMBOL	CRIP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccaagtgtccaagtgcacaaggaggtgtacttcgccgagaggggtgacctctctgggcaaggactggcatcggccctgcctgaagtgcgagaaatgtgggaagacgtgacctctggggccacgtgagcacgaaggcaaacctactgcaaccaccctgtacgcagcc
atgtttgggcctaaaggctttggcgggcgagccgagagccacacttcaagtaa