

C11orf67-chromosome 11 open reading frame 67 Gene

Reference: PTXBC002752

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0366 protein C11orf67; C11orf67; CK067; PTD015; mth938 domain-containing protein; adipogenesis associated Mth938 domain-containing protein; adipogenesis associated Mth938 domain containing		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 67
GENE ID	28971
ORIGIN SPECIES	Human
ACCESSIONS	BC002752
VECTOR	pDONR223
NCBI SYMBOL	C11orf67
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactccctgaaattgcttctatcatggggcgaatgaaagtaaaaggctctatacaacctataaggactgcaaagtatggccagggggtagtcggacttgggattggagagaacaggaaactgagcattctcctggtgtgcagcctgcagatgtgaaggaagttgtgagaaggg
tgtacagactcttgaattggccgaggatgagtgaggccttgaaggtgccttcatcaactgtggagtacctcaagaacatggcattgatgtgcgggtcctccagacagagcaggcagtgaggagtataatgccttgggtgcccaaggggtcagggtgggaggtgtcttcattccacct
gctga