

## IFNAR2-interferon (alpha, beta and omega) receptor 2 Gene

Reference: PTXBC002793

|                |                                                                                                                                                                                                                                                                |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                            | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | IFN-R; IFN-alpha-REC; IFNABR; IFNARB; IMD45; interferon alpha/beta receptor 2; IFN-R-2; IFN-alpha/beta receptor 2; human interferon alpha/beta receptor; interferon (alpha, beta and omega) receptor 2; interferon alpha binding protein; interferon receptor; |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | interferon (alpha, beta and omega) receptor 2                                                     |
| <b>GENE ID</b>          | 3455                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC002793                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | IFNAR2                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgctttgagccagaatgcctcatcgtcagatcacttaattgttctcatggtgtatcagcctcgtgtttggtattcatatgattgcctgattacacagatgaatctgcacttcaagatatcattgcgaaattccgggtccatctatcatgggaattaaaaaacactccattgtaccaactca  
ctatacattgctgtatacaatcatgagtaaacagaaagatttgaagggtgttaagaactgtgcaaataccacaagatcattttgtgacctcacagatgagtggaagcacacagagcctatgtcaccgtcctagaaggattcagcggaacacacagttgttcagttgctcacacaatttct  
ggctggccatagacatgtctttgaaccaccagagttgagattgttttaccacacacattatgtgatggtgaaattccatctattgttgaggagaattacagtttgatttatctctcgtcattgaagaacagtcagaggggaattgttaagaagcataaacccgaataaaaaggaaacatga  
gtggaaattcaccttatcattgacaagttaattccaacacgaactactgtgtatctgtttatttagagcacagtgatgagcaagcagtaataaagtctcccttaaaatgcacctccttcacctggccaggaatcagaatcagcagaatctgccaaaataggaggaataattactgtgttttga  
tagcattggtcttgacaagcaccatagtgacactgaaatggattgttatatatgcttaagaatagcctcccaagcttgaggcaaggctctcactaagggtggaatgcagtgcttattcacaggtgcagtcataatgcactacagtcgtgaaactctgagctcaaacagtcgtctgcctaa  
gcttccccagtagctgggattacaagcgtgcacccctgtgccccagtgattaa