



Reference: PTXBC002801

|                |                                                                                                                                                                         |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                     | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                    | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | P47; UBX1; UBXD10; UBXN2C; dJ776F14.1; NSFL1 cofactor p47; NSFL1 (p97) cofactor (p47); SHP1 homolog; UBX domain-containing protein 2C; p97 cofactor p47; NSFL1 cofactor |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | NSFL1 (p97) cofactor (p47)                                                                       |
| <b>GENE ID</b>          | 55968                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC002801                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | NSFL1C                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

atggcggcgcgagcgcacaggaggcgctgaggggatctgtgcggcgccggcgccgaggaggaccgggcccgcctcttctcgaagtcggccggcttggaacttgcagatcgcgtacgcagcctttatgaggacggagggaatgaagacatttgtaccatttcgcaggcaaccccag  
ttcattgtctcagacggccagcccccagtgataatgagctgacatcttcagacagctcattatgcacgaagatgaagatgagggagggaagaggagccagagctttatgtcggcgccgcagagagaaggaggacagatttgttcggccccaagagaagaatgccaacgcagctg  
ctggatgattctctaaaggctgcgaagagcgcagctgagctgctgtggagcgcagtgacgaagccggcgagcagcagctaaacccagacatttgcaggagctgtctaccgccttggggagcaccagagagagagctgctgattgttcgcagagaagaaagcagcattccac  
caagatgtctgattgattgaacctctggagagctgattacgcctctgataatgggaacttcagagctaccacagccatccatgcacgaatttgcagctgtatccgcagcggggagctgcagcagagctctggagctgctacccgtgttcagctgagctgattgattgagacac  
attcgagcaggagactttgtgcagcccaagcagctctcaaaagctcagctggcgaggctgcagaagctggcgacatccccccagctgtttgatattacagctctccagcccaacagcagaanaatgaagccaaagcagctctcattctatcagatcagcaatccactcaacaac  
atcaaaatctgcgttcagacgcggggagcgtctgctgcagaatfttaaccaacccagagatcagcagcatcgactctcattcgtggatcccgccggccagcattgcgtccacagctttatctctatgactattctccgaacaagaagctggtctgatgagagccagacctgaaggaa  
gccaacctctcaatcgtctatctctcagcgtcttaacataa