

SUPT4H1-suppressor of Ty 4 homolog 1 (*S. cerevisiae*) Gene

Reference: PTXBC002802

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPT4; SPT4H; SUPT4H; Supt4a; transcription elongation factor SPT4; DRB sensitivity-inducing factor 14 kDa subunit; DSIF p14; small hSpt4 subunit; suppressor of Ty 4 homolog 1; SPT4 homolog, DSIF elongation factor subunit		

For research use only.

GENE DESCRIPTION	suppressor of Ty 4 homolog 1 (<i>S. cerevisiae</i>)
GENE ID	6827
ORIGIN SPECIES	Human
ACCESSIONS	BC002802
VECTOR	pDONR223
NCBI SYMBOL	SUPT4H1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctggagacgggtgccgaaggacctgcggcatctgcgggctgtttgtgttcgtggtcaagactatagaccagtttgaatatgatggttgacaattgtgatgcatactacaaatgaagggttaaccgagagatggtatatgactgcactagctcttctttgatggaatcattgcg
atgatgagtcagaggacagctgggtctccaagtggcagcagtcagtaactttaagccaggtgtatatafctggttcagtcactggtgcctgccccaaggaaatcgtcgaggagctgaaaagtcgaggagtgccctacaaatccagagacacagctataaagacctag