



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GJB1-gap junction protein, beta 1, 32kDa Gene

Reference: PTXBC002805

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMTX; CMTX1; CX32; gap junction beta-1 protein; GAP junction 28 kDa liver protein; connexin-32; gap junction protein, beta 1, 32kDa; gap junction protein beta 1		

For research use only.

GENE DESCRIPTION	gap junction protein, beta 1, 32kDa
GENE ID	2705
ORIGIN SPECIES	Human
ACCESSIONS	BC002805
VECTOR	pDONR223
NCBI SYMBOL	GJB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaatggagacggtttgataccctgtcattggcgtgaaccgcgcatctactgcgcattggccgagatgtgctctgcgtcatctctcatcagaaatcagtgctgctggtggtgctgcagagagtgctgggggtatgagaaatctctctcatctgcaacacactcagcctggctgcaacag
cgtttgtatgaccaattctcccacatcccatgtgcggctgtgggtccctgcagctactcctgattccaccccgctctctctgtggccatgcacgtggtcaccagcaacacatagagaagaaaaatgctacggctgaggggccatgggggacccctcacacctggaggagggtgaagaggc
acaggtgccacatctcaggagacatctgggtgacatctgcatcagcgtggttcgcgctgtgtttagagccgctctcatgtatgctcttttactgctcactacctgcctatgcatggctggcctgtcagatgcgacgtctaccctgccccaacacatgtgactgtctgtgtccgccccac
cgagaaaaacgtctaccggtctcatgctctgctctgcctctgcatcatctcgaatcattgagcggcggctggtgttactcatcactccggcgctgtgccccgcgagcccgacccctcccaatccaccttccgccgaagggtcgggtctggccaccgctctcactgaataacagcaga
atgatatcaacaagctctcgaatgacacagatgctccctcgaagacatactctcccgacccctgcacccggcgtggcgtgctgaaagagagcagccctgcgcctgcgcctctga