

MS4A1-membrane-spanning 4-domains, subfamily A, member 1 Gene

Reference: PTXBC002807

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Bp35; CVID5; LEU-16; MS4A2; B-lymphocyte antigen CD20; B-lymphocyte cell-surface antigen B1; CD20 antigen; CD20 receptor; leukocyte surface antigen Leu-16; membrane-spanning 4-domains, subfamily A, member 1; membrane spanning 4-domains A1		

For research use only.

GENE DESCRIPTION	membrane-spanning 4-domains, subfamily A, member 1
GENE ID	931
ORIGIN SPECIES	Human
ACCESSIONS	BC002807
VECTOR	pDONR223
NCBI SYMBOL	MS4A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaacaccagaaatcagtaaatgggactttcccgagagccaatgaaaggccctattgctatgcaatctggtccaaaccactctcaggaggatgtcttactggtgggccccacgcaaagcttctcatgagggaatctaagactttgggggctgtccagattatgaatgggctct
tccacattgccctgggggtcttctgatgatccagcaggatctatgcacccatctgtgtgactgtgtgtgacctctctggggaggcattatgtatatttccggtacactcctggcagcaacggagaaaaactccaggaagtgtttgtcaaaaggaaaaatgataatgaattcattgagcc
tctttgtgccatttctggaatgatttctcaatcatggacatacttaatatataaattcccattttttaaaatggagagctgaattttattagagctcacacacatatataacatatacaactgtgaaccagctaatccctctgagaaaaactccccatctaccaatactgttacagcatacatctc
tgttcttgggcattttgtcagtgatgctgcttcttccaggaaacttgaatagctggcatcgttgagaatgaatgaaaagaacgtgtcagacccaatctaatactgtctcgtcagcagaagaaaaaagaacagactattgaataaagaagaagtgttgggctaactga
aacatcttccaaccaagaatgaagaacattgaaattattccaatccaagaaggaagaagaagaacagagacgaacttccagaacctcccaagatcagggaatctctaccaatagaataatgacagctctcttaa