

SSX2-synovial sarcoma, X breakpoint 2 Gene

Reference: PTXBC002818

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein SSX2; CT5.2; CT5.2A; HD21; HOM-MEL-40; SSX; cancer/testis antigen 5.2; cancer/testis antigen family 5, member 2a; sarcoma, synovial, X-chromosome-related 2; synovial sarcoma, X breakpoint 2; synovial sarcoma, X breakpoint 2B; tumor antigen HOM-MEL		

For research use only.

GENE DESCRIPTION	synovial sarcoma, X breakpoint 2
GENE ID	6757
ORIGIN SPECIES	Human
ACCESSIONS	BC002818
VECTOR	pDONR223
NCBI SYMBOL	SSX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacggagacgacgctttgcaaggagaccacgggttggtgctcaataaccagagaagatccaaaaggccttcgatgatatgccaaatactctctaaggaaagatgggaaaagatgaaagcctcggagaaaatctctatgtgtatatgaagagaaagtatgaggctatgactaaact
aggtttcaaggccaccctccacctttcatgtgtataaacgggcccgaagactccagggaatgattggataatgacctaacctgggaatcaggttgaacgtctcagatgactttcggcaggctccagggaatctcccgaagatcatgcccaagaagccagcagagggaagaaat
gattcggaggaagtgcagagcatctggcccaaaatgatgggaagagctgtgccccgggaaaaccaactactctgagaagattcacgagagatctggaataggaggcccaagaaaagggaagagagacgcggaacagctcatcggtggagcagtcagaacacacaca
acattggtcgattcagtttgcacttctatgggtgcagttcatgtacccccaaacaattacacacacagggacccaaaagggggaacatgcctggaccacagactgcgtgagagaaaacagctgtga