

TMEM11-transmembrane protein 11 Gene

Reference: PTXBC002819

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf35; PM1; PMI; transmembrane protein 11, mitochondrial; transmembrane protein 11		

For research use only.

GENE DESCRIPTION	transmembrane protein 11
GENE ID	8834
ORIGIN SPECIES	Human
ACCESSIONS	BC002819
VECTOR	pDONR223
NCBI SYMBOL	TMEM11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcttggggaaggagcgtcttggccggcagcagtgggcggcagcggcggagagaggggtgagcttgcggccacagactgctacattgtgcatgagatctacaatggggagaatgcccaagaccagtttgagtacgagctggagcaggccctggaagcccgatcaagta
cattgtgattgagccactcgcattggcgacgagacagcccgctggatcacctgggcaactgcctgcacaagacggcgtgctggcgggcaccgctgctcttcacccggtggcgtgccccttagattatccactacattccctgcccgtggtgtgctgagcctggcctgctgca
ccctctatgggatctcctggcagtttgacccttgcctgcaagtaacagtgagtgacgacgacctataaactgtcgcgcctgctcgcacacactcactcctccaccccggtggtgctggtccggaaggacgacctgcacagaaagagactgcacaacacgatagcactggccgcccgtggt
gtactgtgtaagaagatttacgaactctatgccgtatga