

MGC4294-hypothetical protein MGC4294 Gene

Reference: PTXBC002831

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized MGC4294		

For research use only.

GENE DESCRIPTION	hypothetical protein MGC4294
GENE ID	79160
ORIGIN SPECIES	Human
ACCESSIONS	BC002831
VECTOR	pDONR223
NCBI SYMBOL	MGC4294
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgtctctggcccaggagtggtgggggtgcctcacacctgctgtgtggggcagagtctctctgggtcggggcctcgtcggcgatctgtccatgtcgtgtgatgtcgtcagccccgcgtcaccgagagcctcttatgggccaccaagtccacagccactggccatggggca
ccctgcagggcctctggacacggagacctgtgttctcagccagtggtgattttgtgtccagggaacacctggccaagtctgggttcaattttgtgtcacagttaa