

C7orf23-chromosome 7 open reading frame 23 Gene

Reference: PTXBC002837

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein C7orf23; C7orf23; MM-TRAG; MMTRAG; transmembrane protein 243; MDR1 and mitochondrial taxol resistance associated; MDR1- and mitochondrial taxol resistance-associated protein; transmembrane protein 243, mitochondrial		

For research use only.

GENE DESCRIPTION	chromosome 7 open reading frame 23
GENE ID	79161
ORIGIN SPECIES	Human
ACCESSIONS	BC002837
VECTOR	pDONR223
NCBI SYMBOL	C7orf23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggacttgcaccaggacacacggcaccagtgccctggacaacagacctctgttggggagagctccgccaaggatcgaatcatcaatttagttgtggcagcttaacatcctattgattctagtaacgctgataagtgtttgtttccctcaactacacccaaaccgttgatataatc
tttgcgtcgtcatctctttgagtagtattactgctgcatactatctactggtatcgacaaggagacttagaaccgaaattagaagctaatttactatcatatatttctatcatcatgtttgtgtatgtgcaaacctgtacttccatgatgtgggaaggtga