

C1orf35-chromosome 1 open reading frame 35 Gene

Reference: PTXBC002843

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MMTAG2; multiple myeloma tumor-associated protein 2; multiple myeloma transforming 2; chromosome 1 open reading frame 35		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 35
GENE ID	79169
ORIGIN SPECIES	Human
ACCESSIONS	BC002843
VECTOR	pDONR223
NCBI SYMBOL	C1orf35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcggtccagtcgtggagcgtgcgcggcggcaggaccagttcaactgggagcgtgaagactgacaagcagcgggagaactacgtggcaactcgtgatggcgcggtaggcccgtggcagaaggccgcgacctcacctggtacgccaaggccggcgccatg
cgccggcccgagccgcgagaggaactggcagccgtgcgggagggcgagcgcgagcgtctgctggccgcccttggttacaagaactgaagaagcagcccccggcctgagcaaggaggacttcgcggaggtctgcaagcgggaaggagcgacccccgaggagaaggc
gtggaccgctgctgggctgggagcgaagtggctccgtggcgcgtggcgtatgtccgagaggacaaggagccgccaactgggctgtctgtttcacgcatcaccgcgtagagagcggcggccgggacctggcagcctcgccaggaggagccggcgccg
aggatcagacggaaagcagttgtgagagccacaggaagcaaggaaggaagaaaaagaaaaagaggaacacaaagaaagagaagaagaagacaaagagcagcggccagctgagccacctctctcccatctctgagaggccagccaccacc
atgactccgactccaactccccctgctgaagaggaggagcggggacacagtggggacaggaggagcccgtctcgcaggtggcatgacagaggctctgagggcctga