

C2orf28-chromosome 2 open reading frame 28 Gene

Reference: PTXBC002846

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C2orf28; APR--3; APR-3; APR3; HSPC013; PRO240; p18; all-trans retinoic acid-induced differentiation factor; apoptosis related protein APR-3; apoptosis-related protein 3; all-trans retinoic acid induced differentiation factor		

For research use only.

GENE DESCRIPTION	chromosome 2 open reading frame 28
GENE ID	51374
ORIGIN SPECIES	Human
ACCESSIONS	BC002846
VECTOR	pDONR223
NCBI SYMBOL	C2orf28
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgcatgccgttgctgcctgaatcagaagggcaccatctggggctggatccagaactgttctctggagaccctgtccaaacttcatcaggcacataccactgtcatcatagacctgcaagcaaacccctcaaaggtgacttgccaacaccttcctggtttactcagctcc
agactctgatactgccacaacatgtcaactgtcctggaggaattaatgcttggaatactatcacctcttatatagacaacaaatctgtcaaggcgaagaacatttgcataaacactggggaccagaaaatgtctctgagaatggatcttgtgtacctgatggtccaggtctttgcagtgtgt
ttgtgctgatggtttccatggatacaagtgatgcgccagggctcgttctcactgcttatgttcttcgggacttgggagccaccactctatccgtctcattctgcttgggcgaccagcgccgaaaagccaagacttcatga