

C3orf60-chromosome 3 open reading frame 60 Gene

Reference: PTXBC002873

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf60; 2P1; E3-3; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 3; NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, assembly factor 3; NADH dehydrogenase (ubiquinone) complex I, assembly factor 3; nuclear protein E3-3; NADH:ubiqu		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 60
GENE ID	25915
ORIGIN SPECIES	Human
ACCESSIONS	BC002873
VECTOR	pDONR223
NCBI SYMBOL	C3orf60
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaccgctctcgctacgttagctgtaccgagcgcgacctcgtcgtcgtgtccgcccgtgagcttccctgggccccgcggcgaggcgatcgctctcggcgcgatgacgagctgtatcagcggacgcgcacatctctgtgcaacgcgagccgctcaggcaatgtacat
cgacagctacaacagccgcttcatgataaacggaaaccgctgctcggccctgctcgtctcgtccgcactcgggtggtcagtggaacgtgggatccaccaggacatcaccgaagacagcttttccctctctgtgtgagccccgagatagagatcgtgtgtgtgggactg
gagaccggaccgagagcgtcagtcctcaggtgctcaagccatgagcagcggggcattgctgtggaatgaggacacgcccattgctgtccacctcaacttctgtgtcatgaagccgagtaactggagctgctctcatccctccaccaggaggacttcacttactcttgg
gccaagctgctcaatga