

SETD4-SET domain containing 4 Gene

Reference: PTXBC002898

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf18; C21orf27; SET domain-containing protein 4; SET domain containing 4		
For research use only.			

GENE DESCRIPTION	SET domain containing 4
GENE ID	54093
ORIGIN SPECIES	Human
ACCESSIONS	BC002898
VECTOR	pDONR223
NCBI SYMBOL	SETD4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagaaggaaaggagaaacagccgatcagaagacgaaactctgcggaggttctgaatcaaggagtgatgagagccacaagctgaattatagagctgagggaagtggctgaaagctaggaagttcaagattcaacttagcgctgtgtttccaggtacaggaaga
gggctgatgagtaaacatccctgcagaggggacagatgatttctgttgcctgagagttgcctgctcaccacggacacagtgattcgaagctacttaggggcatacattactaagtggaagcctcctccatctctctgtggcgctgtgcaccttttagttcagaaaagcatgctgggca
ccgatctcttggaaagccttacctggagattttacccaaggcgtatacctgcccgtttgttggagccggagtggtgaaccttctcccaatctttaaagcaaaagctgaagagcagagagccacgtgcaggagttcttcttctccagagacttttctctctgcagcctctgtttgc
ggaggtgttgacagcatctcagctacagtgccctgctgtgggttgggtgcacgtcaacaccagagccgtgtacgtgagggccaggcagcgggaatgccttctgcagagccggacacctgtgcactcgtcctgacgtggacctgtgaatcatagcccacatgtccaggtaaaagc
agcgtttaatgaagaactcattcttacgaaattagaacgacttcacgttggagaaagcatgaagaggtattcattgttacggccctcacgataatcaacggctgttctggaatacggattgttctgtccataatcctcatgtgtttatgtctcaagaggttggaatcaactttgttctaa