

MAD2L1BP-MAD2L1 binding protein Gene

Reference: PTXBC002904

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAD2L1-binding protein; caught by MAD2 protein; MAD2L1 binding protein		
For research use only.			

GENE DESCRIPTION	MAD2L1 binding protein
GENE ID	9587
ORIGIN SPECIES	Human
ACCESSIONS	BC002904
VECTOR	pDONR223
NCBI SYMBOL	MAD2L1BP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgccggagggtgtgtctcagcgcagtcctgattggagtggtatgagaagtcgaagaactcagcctccagatagaactctgagacaagctctacgaggaacctctcaacgttcggaggcctttgcccaagagactgcatggtaccagtgtgtttc
ctgggcctgtgagccagggaaggctgtgtcagtttactgtgaacttctaaagcatatcatgtatcaacgccagcagctccctctgcccataaacagcttaagcacttttaccgaaaacctctcccaggcagagggagatgctgaagaagaacctgggcccactgaggtgagcagc
aggaaatgccaacaagccctggcagaactggagagtgctcagccacctggaggacttctttgacggacactagtagccgagtgctgattctcttggggcgaatgccctaagcccaaggagttctatgaactcgaactgtctctgtggccccctacagcgtggaccagagcctga
gcacagcagctgtttgcgcgctctctccgagccatattcagctgagccttttagcagcttcaggtctctccactatgggcaccgtctgcatggcacaggacaccgcaactgtggagaagattgtttcgacccaagctcaactatcagtgcccagccggggccataaactgactg
tgacctgtcatgtggcagaccttccatccgaaccacgcttgggaagactacattgtgtccagccaccagtgacatttaaggcttcgcgagtgga