

RPRM-reprimo, TP53 dependent G2 arrest mediator candidate Gene

Reference: PTXBC002908

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	REPRIMO; protein reprimo; candidate mediator of the p53 dependent G2 arrest; reprimo, TP53 dependant G2 arrest mediator candidate; reprimo, TP53 dependent G2 arrest mediator candidate		

For research use only.

GENE DESCRIPTION	reprimo, TP53 dependent G2 arrest mediator candidate
GENE ID	56475
ORIGIN SPECIES	Human
ACCESSIONS	BC002908
VECTOR	pDONR223
NCBI SYMBOL	RPRM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatccggccctaggcaaccagacggacgtggcgccgtgttctctggccaacagcagcagggcgctggagcagccgtgcgctgctgcacccaggcgctccgtggtagccgacgacggcttcggaggaggagggccggagcagcgtagcctgtacataatgcgcgtggtgcag
atcgcggtcatgtgcgtgctctcactcaccgtgggtcttggcgtcttctctcggctgcaatctgctcatcaagtccgagggcatgatcaacttctctgtaaggaccggaggccgtctaaggaggtggagcggtgtgctgtgggcctactga