

VKORC1-vitamin K epoxide reductase complex, subunit 1 Gene

Reference: PTXBC002911

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EDTP308; MST134; MST576; VKCFD2; VKOR; vitamin K epoxide reductase complex subunit 1; phyloquinone epoxide reductase; vitamin K dependent clotting factors deficiency 2; vitamin K1 2,3-epoxide reductase subunit 1; vitamin K1 epoxide reductase (warfarin-se)		

For research use only.

GENE DESCRIPTION	vitamin K epoxide reductase complex, subunit 1
GENE ID	79001
ORIGIN SPECIES	Human
ACCESSIONS	BC002911
VECTOR	pDONR223
NCBI SYMBOL	VKORC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcagcactgggggagccctggctgggtgcgctcctttgcctgacgggcttagtgcctcgcctacgcgctgcacgtgaaggcggcgcgccgggacccgggattaccgcgcgtctgcgacgtgggcaccgccatcagctgttcgcgcgtctctccaggtggg
gcaggggtttcgggctggtggagcatgtgctgggacaggacagcatcctcaatcaatccaacagcatattcgggtgcatcttctacacactacagctattgttaggttgcctgcggacacgctgggcctctgtcctgatgctgctgagctccctggtgtctctcgtggttctgtctacctggcct
ggatcctgttctcgtgctctatgatttctgcatgtttgtatcaccacatgctatcaacgtgagcctgatgtggctcagtttccggaaggtccaagaacccagggaaggcctaagaggcactga