



Reference: PTXBC002914

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRPL-2; WAS2; WASPIP; WAS/WASL-interacting protein family member 1; WASP-interacting protein; Wiskott-Aldrich syndrome protein interacting protein; protein PRPL-2; testicular tissue protein Li 226; WAS/WASL interacting protein family member 1		

For research use only.

GENE DESCRIPTION	WAS/WASL interacting protein family, member 1
GENE ID	7456
ORIGIN SPECIES	Human
ACCESSIONS	BC002914
VECTOR	pDONR223
NCBI SYMBOL	WIPF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]