



Reference: PTXBC002917

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPR172A; BVVLS2; D15Erd747e; GPCR41; PAR1; RFT3; RFVT2; hRFT3; solute carrier family 52, riboflavin transporter, member 2; G protein-coupled receptor 172A; PERV-A receptor 1; porcine endogenous retrovirus A receptor 1; riboflavin transporter 3; solute ca		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 172A
GENE ID	79581
ORIGIN SPECIES	Human
ACCESSIONS	BC002917
VECTOR	pDONR223
NCBI SYMBOL	GPR172A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]