

TSPAN14-tetraspanin 14 Gene

Reference: PTXBC002920

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DC-TM4F2; TM4SF14; tetraspanin-14; tetraspanin similar to TM4SF9; transmembrane 4 superfamily member 14; tspan-14; tetraspanin 14		

For research use only.

GENE DESCRIPTION	tetraspanin 14
GENE ID	81619
ORIGIN SPECIES	Human
ACCESSIONS	BC002920
VECTOR	pDONR223
NCBI SYMBOL	TSPAN14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcactattatagatacttaacccaaggtcagctgctgtgtaacaagtaacctcttttcagctacaacatcatcttctgggtgtgtgtctcgacctcaccagaagtgacccggatgcatggaatcgaccctgtggtgctgtctgatgtgtggcggtgtgtgtgttcaccctgggttcgcccgtgctgtggggctctgctggggagaatatctgtgtctcaactttttctgtggcaccatcgtgtctcatcttctctgagctggctgtgtggcgtgtgtgtgctgtgcttctgttccaggactgggtgagggaccggttcggggagttcttcgagagcaacatcaagtcctaccgggacgatactgatctgcaaaacctcatcgactccttcagaaagctaaccagtgctgtggcgcatatggcctgaagactgggacctcaacgtctacttcaatgtcagcgggtgccagctacagccgagagaagtgccgggtcccttctctgctgtgtgccagatcctgcgcaaaaagtgtgaacacacagtggtgataatgatgtcaggatcagctgaagagcaagtggtgatgagtccttcacgaaaggctgcatccaggcgtggaagctggctcccgcggaacatttacattgtggctggcgttcatcgccatctgctgttgagatattggcatcttctggcaaggacgctgatctcagacatcgagcgagtgaaggccggccatcactctga