



Reference: PTXBC002922

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	polycomb complex protein RING1; E3 ubiquitin-protein ligase RING1; RING1A; RNF1; really interesting new gene 1 protein; ring finger protein 1		

For research use only.

GENE DESCRIPTION	ring finger protein 1
GENE ID	6015
ORIGIN SPECIES	Human
ACCESSIONS	BC002922
VECTOR	pDONR223
NCBI SYMBOL	RING1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr