



Reference: PTXBC002946

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	miro domain-containing protein C1orf89; C1orf89; REM2- and Rab-like small GTPase 1; Rem/Rab-Similar GTPase 1; REM2 and RAB like small GTPase 1		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 89
GENE ID	79363
ORIGIN SPECIES	Human
ACCESSIONS	BC002946
VECTOR	pDONR223
NCBI SYMBOL	C1orf89
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atgccagacacccctgctccggctgtgctggtgtcccaacgtgacacgagatgccgaggccaaagatgactcgtctcatctctgcgaagaaccgccggcggtgttgggctgctgagcggccagctgctctgcgcctgtgtccatgacactgccagctacaagaactttgtgtccgggaagagtggtgtgggccaagacggcgctggtgccaaagtctggtccctggaggtgctgtggtgccaccacgggacaccgggcatccagaccaccgtggtatfttggccagccaagtctgcaaggccagcagccgtgtctgcatgttttgagtctgggagctggagagctgcactcaaaaatctgatcatatgctgctgctgcatggagaacacagatgctctctctctctctctctcaatgaccgtgctcctcttggaaagacacctcctggacagctggcccgcatagacaggtgagggccctgggtgtctgagatggtcaatggctcccaattggacagtacatgcacaggacgtgcccggagcgggacatccctggcaggcgctgggagctgcccctgctacgggtgaagagtgctggccggcgccggtggtgatggccgcacatgcccggcgccggcggtggtggcgactgtgccatactcaatggccttctgtagcagctgtggcacaggacacagctggcctgctgctcctcccaaccccccaagagatgctctgaatga

Page 1/1