

PSMA6-proteasome (prosome, macropain) subunit, alpha type, 6 Gene

Reference: PTXBC002979

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IOTA; PROS27; p27K; proteasome subunit alpha type-6; 27 kDa prosomal protein; PROS-27; macropain iota chain; macropain subunit iota; multicatalytic endopeptidase complex iota chain; prosomal P27K protein; proteasome (prosome, macropain) subunit, alpha typ		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) subunit, alpha type, 6
GENE ID	5687
ORIGIN SPECIES	Human
ACCESSIONS	BC002979
VECTOR	pDONR223
NCBI SYMBOL	PSMA6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcccggtgtccagcgccggttttgaccgccacattaccatttttaccgccggtcggctctaccaagtagaataatctttaaaggctattaaccagggtggccttacatcagtagctgtcagaggggaaagactgtgcagtaattgtcacagaagaagacacgtgacaaattattggatt
ccagcacagtgactcacttattcaagataactgaaacattggtgtgtgatgacccgaatgacagctgacagcagatccaggtacagagggcacgctatgaggcagctaacctggaatacaagatggctatgagattcctgtggacatgctgtgtaaaagaattgccgatatctcaggt
ctacacacagaatgctgaaatgaggcctctgtgtgtgtatgatttaattgggtatagatgaagagcaaggccctcaggtatataagtgatcctgcaggttactactgtgggtttaagccactgcagcgggagttaaacaaactgagtcacacagcttcttgaaaaaaagtgaagaaga
aatttgattggacatttgaacagacagtggaactgcaattacatgcctgtctactgttctcaattgattcaaaccttcagaataagaagttggagtagtgacagttgaaatcctaattcaggattcttacagaagcagagattgatgctcacttgttgccttagcagagagagactaa