

ZNF576-zinc finger protein 576 Gene

Reference: PTXBC002981

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 576		
For research use only.			

GENE DESCRIPTION	zinc finger protein 576
GENE ID	79177
ORIGIN SPECIES	Human
ACCESSIONS	BC002981
VECTOR	pDONR223
NCBI SYMBOL	ZNF576
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggacccgaacctgaagagaacatgaagcagcaggtaccacgaagagagaagtcacccagagccagaggcgaacatctgccacctggggggcccggaagtgcaccgctgcctcatcaccttcagattccaagttccaggagcgtcacatgaagcgggagcaccag
cggacttcgtggccagaagctgcaggggtctctctcatctgttcacctgcgcgcctctctctccaaagccctaatcaccaccagcgcagccacggtccagccgccaagcccacctgcgggttgcaaccactactgccagcccaccttccctgtctgactgtggcaag
acctttggcgaggctgtttctctgagggcgccagccagatgcatgaggtccgtgcccctctggcaccttcgctgcacagatgcgggtcaggactttgtcaggaagcagggctgcatcaactacattcggcatgcccggggggagctctga