

PAFAH1B3-platelet-activating factor acetylhydrolase, isoform Ib, gamma subunit 29kDa Gene

Reference: PTXBC003016

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PAFAHG; platelet-activating factor acetylhydrolase IB subunit gamma; PAF acetylhydrolase 29 kDa subunit; PAF-AH 29 kDa subunit; PAF-AH subunit gamma; PAF-AH1b alpha 1 subunit; PAFAH subunit gamma; platelet-activating factor acetylhydrolase 1b, catalytic s		

For research use only.

GENE DESCRIPTION	platelet-activating factor acetylhydrolase, isoform Ib, gamma subunit 29kDa
GENE ID	5050
ORIGIN SPECIES	Human
ACCESSIONS	BC003016
VECTOR	pDONR223
NCBI SYMBOL	PAFAH1B3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtggaggagaaacccagccagcaagccacgcccgtgacggacgtacagggcgacggcgctggatgtccctgcaccatcggtctggtctgacagcaagataaggaacccgaagtcgtcttcacggtgggactcctgtggtccagctcatgcaccagtgcgagatctggcg
cgagctctctctctctgcatgcaacttaacttggcattggtgacggcacacagcatgtactgtggcggtggagaatggggagctggaacacatccggcccaagattgtgtgtgtctggtgggacccaacaaccacggacacacagcagacagtgactggtggcatcaagg
ccattgtgcaactggtgaatgagcgacagcccccagcccggttgggtgctggcgctgtccgcgagggccaacatcccaacccacttgggagaagaacggacaggtgaacgagctgtacggcgcgacgtggctggccaccctcgggccacttcttagatccgacccctggct
ttgtgactcagatggcaccatcagccatcatgacatgtatgattacctgcatctgagccgctgggtacacacctgttggcggtctgcactccctgcttctgctgtctgtggccaagaccagggccaaggtgctccctgctggagcccgcacccaa