



Reference: PTXBC003060

|                |                                                                                                                                |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                            | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | dJ20N2.5; plasma alpha-L-fucosidase; alpha-L-fucosidase 2; alpha-L-fucoside<br>fucohydrolase 2; fucosidase, alpha-L- 2, plasma |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | fucosidase, alpha-L- 2, plasma                                                                   |
| <b>GENE ID</b>          | 2519                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC003060                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | FUCA2                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

[illegible]