

ILF3-interleukin enhancer binding factor 3, 90kDa Gene

Reference: PTXBC003086

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CBTF; DRBF; DRBP76; MMP4; MPHOSPH4; MPP4; NF-AT-90; NF110; NF110b; NF90; NF90a; NF90b; NFAR; NFAR-1; NFAR2; TCP110; TCP80; interleukin enhancer-binding factor 3; M-phase phosphoprotein 4; double-stranded RNA-binding protein, 76 kD; dsRNA binding protein N		

For research use only.

GENE DESCRIPTION	interleukin enhancer binding factor 3, 90kDa
GENE ID	3609
ORIGIN SPECIES	Human
ACCESSIONS	BC003086
VECTOR	pDONR223
NCBI SYMBOL	ILF3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgttcatgtctatctctctcctggacacattaagggttgggacaagaagagactcctgtgagatcagcagacctcgtgtttaaggagagaatgcctgcctgcaacactcatggagaggacacttcctggccccacctgctaataacctccaaactctggcaggcgctgcattctc
aacacaggacttccttaacagtttccggtatcagtgaccaccacagttcaggtttctgaatgcagaccaccacaggggatggaftccaccgcccgccctgctgcagatacatgcggcacagcagcagcagtgctggagagaacagctcctgagcagcagccttgagccctca
ccccgtcacgatcctcaccagtgctgatgatccgtcatgctggtgcctgagacagatccagctacaacaggacgacaccgggttcaggaatatgtgctatcacctga