

RRAGD-Ras-related GTP binding D Gene

Reference: PTXBC003088

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAGD; bA11D8.2.1; ras-related GTP-binding protein D; Rag D protein; Ras related GTP binding D		

For research use only.

GENE DESCRIPTION	Ras-related GTP binding D
GENE ID	58528
ORIGIN SPECIES	Human
ACCESSIONS	BC003088
VECTOR	pDONR223
NCBI SYMBOL	RRAGD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagccctggccaggtccacctcacgtgaccagggcctacaaagtgaatactgacatcaacttcgaggtgtttatcataaaggatgggtctgtcagatgaccacaaaattgaaacccaaagagatattcaccagaggcacaacgatgaccttcagatgctggattagaaaaat
tcacctcagcttttatctgacaagcatatgatcattcaatattgaaagcttttagcaaaagtgttcagaaactgattccacaactcccaactctggagaatttgcgaacatctttatctcaaatctggaattgaaaaggcatttctattgatggtcagtaaaatttatattgcaactgatagtactcc
gggtgatatgcaaacctatgagctctgctgatatgatatggtttatgacatctctgtatttgggtctcaaaagaagatggagcaggaacccctatgacaaggaaatccacagccatcataaagcttaataatacaaccgtgctttattttaaagaggtgacaaagtctcctgctctgttg
ctttgtcagagaggaagctttgaaagaaaaggcctaattgactataatttcattgctccggaaggccattcatgaagttttgaggtgagaatgaaagtagtaaaatctcgaagggttcagaatcggtgcagaagaaaagagagcccccctaattgggacccctagatgctgctgtag