



Reference: PTXBC003091

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PVRL2; CD112; HVEB; PRR2; PVRR2; nectin-2; herpesvirus entry protein B; poliovirus receptor-like 2; poliovirus receptor-related 2 (herpesvirus entry mediator B); nectin cell adhesion molecule 2		

For research use only.

GENE DESCRIPTION	poliovirus receptor-related 2 (herpesvirus entry mediator B)
GENE ID	5819
ORIGIN SPECIES	Human
ACCESSIONS	BC003091
VECTOR	pDONR223
NCBI SYMBOL	PVRL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]