

TBCD-tubulin folding cofactor D Gene

Reference: PTXBC003094

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PEBAT; SSD-1; tfcD; tubulin-specific chaperone D; beta-tubulin cofactor D; tubulin folding cofactor D		

For research use only.

GENE DESCRIPTION	tubulin folding cofactor D
GENE ID	6904
ORIGIN SPECIES	Human
ACCESSIONS	BC003094
VECTOR	pDONR223
NCBI SYMBOL	TBCD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr