

ICAM2-intercellular adhesion molecule 2 Gene

Reference: PTXBC003097

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD102; intercellular adhesion molecule 2; ICAM-2		
For research use only.			

GENE DESCRIPTION	intercellular adhesion molecule 2
GENE ID	3384
ORIGIN SPECIES	Human
ACCESSIONS	BC003097
VECTOR	pDONR223
NCBI SYMBOL	ICAM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctcttgcgttacaggaccctgactgtggccctcttcaccctgatctgtctcaggatcggaaggtattcgaggtacacgtgagccaaagaagctggcgggtgagccaaagggtccctcagggtcaactgcagcaccacctgtaaccagcctgaagtgggtgtctgg
agacctcttagataagattctctggacgaacaggctcagtggaacattacttggctcaaacatctccatgacacgggtcccaatgccacttcacctgtccgggaagcaggagtcattgaattccaacgtcagcgtgtaccagcctcgaaggcaggtcatctgacactgcaaccca
cttgggtggtgtggcgaagtccttcaccattgagtgagggtgcccaccgtggagccctggacagcctcaccctctctctgttccgtggcaatgagactctgactatgagaccttcgggaaggcagccctgtcccgaggaggccacagccacattcaacagcagcgtgacagag
aggatggccaccgcaacttctctgctgctggtgctgacttgaatgtctcgggtggcaacatcttcacaacactcagcccccgaagatgttggagatctatgacctgtgtcggacagccagatggtcatcatagtcacgggtgtgtcgggtgtgtgtctgtctgtgacatctgtctg
ctctgttcatcttggccagcacttgcgcagcagcggatgggcacctacgggtgctgagcggcttgaggaggtgccccaggccttcggccatag