

CPSF4-cleavage and polyadenylation specific factor 4, 30kDa Gene

Reference: PTXBC003101

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPSF30; NAR; NEB-1; NEB1; cleavage and polyadenylation specificity factor subunit 4; CPSF 30 kDa subunit; NS1 effector domain-binding protein 1; cleavage and polyadenylation specific factor 4, 30kDa; cleavage and polyadenylation specificity factor 30 kDa		

For research use only.

GENE DESCRIPTION	cleavage and polyadenylation specific factor 4, 30kDa
GENE ID	10898
ORIGIN SPECIES	Human
ACCESSIONS	BC003101
VECTOR	pDONR223
NCBI SYMBOL	CPSF4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggaaatcatcgccagcgtggaccacatcaagttgactggagatcgcggtggagcagcagctggggcgccgctgccttccccggcatggacaagtcggcgctgtgtctgtgaattcttttgaagctgcctgcggcaaaagggggcatgttccgtttcgccacatca
gtggtgagaagacagttgtgtgcaaacactggctgcgtggcctatgcaagaagggggaccagttgtgagttcctgcatgagtagatgacatgaccaagatgccgagtgctacttctactccaagttcggggagtgagcaacaaggaatgtcccttctgcacatcgaccccgagtcgaagat
caaggactgtccttgatgaccgtggccttctgcaagcaggtccctctgcaggcaccggcacacacggagagtcacatctgtgtgaattacctcgtgggattctgcccggagggccctcgtgtaattcatgcaccctcgattgaaactgccatgggaaccacgagcagccccactg
ccgcagcagacagcctccagcaagcagagaaccccgaggtcatcgggtcatgagagtaaacagcagcgcgggcaaccggggaccccgccactggagcaggtcacctgttacaagtgtggcgagaaggacactacgccaacagatgcacaaagggcacttggc
ctttctcagtgacagtga