

NDUFA7-NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 7, 14.5kDa Gene

Reference: PTXBC003102

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B14.5a; CI-B14.5a; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 7; NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 7, 14.5kDa; NADH-ubiquinone oxidoreductase subunit B14.5a; complex I B14.5a subunit; NADH:ubiquinone oxidoreductase subunit		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 7, 14.5kDa
GENE ID	4701
ORIGIN SPECIES	Human
ACCESSIONS	BC003102
VECTOR	pDONR223
NCBI SYMBOL	NDUFA7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtccgccaccggtctcatccagcggctgcggaactggcgtccgggcatgacctgcagggggaagctgcagctacgtaccaggagatctccaagcgaactcagcctctccaagctccctgtgggtctagccacaagctctccaacaattactattgcactcgcgatggccg
ccgggaatctgtgcccccttccatcatcatgtcgtcgcagaaggcgtggtgtcaggcaagccagcagagagctctgtctgtagctgccactgagaagaaggcgggtgactccagctctcccataaagaggtgggagctgtctcggaccagccttacctgtga