



Reference: PTXBC003126

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PCDH-ALPHA2; protocadherin alpha-2; KIAA0345-like 12; protocadherin alpha 2		
For research use only.			

GENE DESCRIPTION	protocadherin alpha 2
GENE ID	56146
ORIGIN SPECIES	Human
ACCESSIONS	BC003126
VECTOR	pDONR223
NCBI SYMBOL	PCDHA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr