

SEPX1-selenoprotein X, 1 Gene

Reference: PTXBC003127

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SEPX1; HSPC270; SELR; SELX; SepR; methionine-R-sulfoxide reductase B1; selenoprotein R; selenoprotein X, 1; methionine sulfoxide reductase B1		

For research use only.

GENE DESCRIPTION	selenoprotein X, 1
GENE ID	51734
ORIGIN SPECIES	Human
ACCESSIONS	BC003127
VECTOR	pDONR223
NCBI SYMBOL	SEPX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgttcgcagcttctcggggcgagggtttccagaatcacttgaacctggcgtttacgtgtgtgccaagtgtggctatgagctgtttccagccgctcgaagtatgcacactcgtctccatggccggcggttcaccgagaccattcagccgcgacagcgtggccaagcgtccggagcaca
atagatctgaagccttgaagggtgtcctgtggcaagtgtggcaatgggttggccacgagttcctgaacgacggcccaagccggcgagtcctcgaattctgaattatcagcagctcgtgaagtgtgtccctaaaggcaagaaacttctgcctccaggggtcactag