

TSPAN15-tetraspanin 15 Gene

Reference: PTXBC003157

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	2700063A19Rik; NET-7; NET7; TM4SF15; tetraspanin-15; tetraspan NET-7; transmembrane 4 superfamily member 15; transmembrane 4 superfamily member tetraspan NET-7; tspan-15; tetraspanin 15		

For research use only.

GENE DESCRIPTION	tetraspanin 15
GENE ID	23555
ORIGIN SPECIES	Human
ACCESSIONS	BC003157
VECTOR	pDONR223
NCBI SYMBOL	TSPAN15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcgcggggactcggagcaggtgcgctactcgcgcgcgttctctacctctggctcaagtttcaactatcatctattccaccgtgttctggctgattggggccctggtcctgtctgtggcctatgcagaggttgagcggcagaaataaaacccctgaaagtgccttctgtgctcca
gccatcatctcatctctcggcgctcgtcatgttcatgtgtctcttcatgtgtgtcgtgcgtccctcgtgacaacctgtaccttccaaagcattcatgtacatcttgggatctgcctcatcatggagctcattggcggtggtggtgaccttccggaaccagaccattgacttctgaac
gacaacattcgaaggaattgagaactactatgatgatctggacttcaaaaacatcatggactttgttcagaaaaagttcaagtgtgtgctggggaggactaccgagattggagcaagaatcagtaaccacagctgcagtgccctggaccctggcctgtgggtgacctacacgtgtg
catcaggaaacacgacagaagttgtcaacacatgtgtgtgctacaaaactatcgacaaggagcgtttcagttgtgcaggatgtcatctacgtgcggggctgcaccaacgccgtgatcatctgtttcatggacaactacacatcatggcggtatcctctgggcatctgttccccagttct
gggggtgctgctgacgtcgtgtacatcacccgggtggagacatcatcatggagcactctgtcactgatgggtccttggggcccggtgccaagccacgcgtggagcggcaggcagcggaatgctgtgtgtacccaattag