



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

POLE3-polymerase (DNA directed), epsilon 3 (p17 subunit) Gene

Reference: PTXBC003166

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CHARAC17; CHRAC17; YBL1; p17; DNA polymerase epsilon subunit 3; CHRAC-17; DNA polymerase II subunit 3; DNA polymerase epsilon p17 subunit; DNA polymerase epsilon subunit p17; arsenic transactivated protein; asTP; chromatin accessibility complex 17 kDa pro		

For research use only.

GENE DESCRIPTION	polymerase (DNA directed), epsilon 3 (p17 subunit)
GENE ID	54107
ORIGIN SPECIES	Human
ACCESSIONS	BC003166
VECTOR	pDONR223
NCBI SYMBOL	POLE3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr