



Reference: PTXBC003168

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ORP10; OSBP9; oxysterol-binding protein-related protein 10; ORP-10; OSBP-related protein 10; oxysterol binding protein like 10		

For research use only.

GENE DESCRIPTION	oxysterol binding protein-like 10
GENE ID	114884
ORIGIN SPECIES	Human
ACCESSIONS	BC003168
VECTOR	pDONR223
NCBI SYMBOL	OSBPL10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]