



Reference: PTXBC003171

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GAA1; hGAA1; glycosylphosphatidylinositol anchor attachment 1 protein; GAA1 protein homolog; GPAA1P anchor attachment protein 1 homolog; GPI anchor attachment protein 1; GPI transamidase subunit; anchor attachment protein 1 (Gaa1p, yeast) homolog; glycop		

For research use only.

GENE DESCRIPTION	glycosylphosphatidylinositol anchor attachment protein 1 homolog (yeast)
GENE ID	8733
ORIGIN SPECIES	Human
ACCESSIONS	BC003171
VECTOR	pDONR223
NCBI SYMBOL	GPAA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]