

MALL-mal, T-cell differentiation protein-like Gene

Reference: PTXBC003179

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BENE; MAL-like protein; mal, T-cell differentiation protein like		
For research use only.			

GENE DESCRIPTION	mal, T-cell differentiation protein-like
GENE ID	7851
ORIGIN SPECIES	Human
ACCESSIONS	BC003179
VECTOR	pDONR223
NCBI SYMBOL	MALL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctgcgcgacccgccaccagctacgcccgtcgacgtgccctcgggggtcgcgctgttcctcacatcccttgccttctctgcgcgagctgatattgggtcttggctcgaccatggtagccgccaccacatagtatacccttgcctgaaggatgggtgatgatg
tctcgtcacctcgtttctcatctccttgatgttcctgtgttcttactgtttggattttacaaaagatttgaatcctggagagttctggacagcctgtaccacgggaccactggcatcctgtacatgagcgtgccgtcctacaagtacatgccagattgtttctgagaaactgctggaccaagaatt
tactacattaattcggcagcctcgtttctgccttcacgccagctgctctacatttccatgccttcagcatctattaccactga