

NGFRAP1-nerve growth factor receptor (TNFRSF16) associated protein 1 Gene

Reference: PTXBC003190

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NGFRAP1; Bex; DXS6984E; HGR74; NADE; protein BEX3; brain-expressed X-linked protein 3; nerve growth factor receptor (TNFRSF16) associated protein 1; ovarian granulosa cell 13.0 kDa protein HGR74; p75NTR-associated cell death executor; brain expressed X-li		

For research use only.

GENE DESCRIPTION	nerve growth factor receptor (TNFRSF16) associated protein 1
GENE ID	27018
ORIGIN SPECIES	Human
ACCESSIONS	BC003190
VECTOR	pDONR223
NCBI SYMBOL	NGFRAP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaaatattcaccaggaaaacgaagagatggagcagcctatgcagaatggagaggagaccgcccttgggaggagggtgaaggccaccagcctgcaggaaatcgacggggacaggtcgcgcgactggccctaatttcgatgggccataccaataggcagatcaatgatggg
atgggtggagatggagatgatatggaaatattcatggaggagatgagagaatcagaagaaaacttagggagctgcagttgagggaattgtctgcgtatccttatggggagctcttaaccacatgaccatcatgatgaatttgccttatgccttga