

TMEM39A-transmembrane protein 39A Gene

Reference: PTXBC003192

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 39A		

For research use only.

GENE DESCRIPTION	transmembrane protein 39A
GENE ID	55254
ORIGIN SPECIES	Human
ACCESSIONS	BC003192
VECTOR	pDONR223
NCBI SYMBOL	TMEM39A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccggtggaaggaggggccctagtcggcaacagctaagccgttcagctttaccttctttgcagactttgtgtggaggctgtggcaatggaacaggcttgagaaacagggaatgtagtgcctattggccttcagtcaccatcacagccttaacacccagggtcctgttcgtcatt
gccaaatctgacttgctgtggatgggagcctactcttgaattcctttttcatctacctgtgtgtgctctttcattcagtaacacattataaaacagtgtgtgtgtatccttacaatcatcctgcttctgtacatcactgaattttcatcctcattgattatcatctggcagcattcatcacagtga
tgcttgcgaggaggctgtatgggctctcatctcagaggtctgtgtgtgtaagctctga