

RPLP1-ribosomal protein, large, P1 Gene

Reference: PTXBC003369

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LP1; RPP1; 60S acidic ribosomal protein P1; acidic ribosomal phosphoprotein P1; ribosomal protein, large, P1; ribosomal protein lateral stalk subunit P1		

For research use only.

GENE DESCRIPTION	ribosomal protein, large, P1
GENE ID	6176
ORIGIN SPECIES	Human
ACCESSIONS	BC003369
VECTOR	pDONR223
NCBI SYMBOL	RPLP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctctgtctccgagctcgctgcatctactcggccctcattctgcagcagcatgaggtgacagtcacggagataagatcaatgccctcattaaagcagccggtgtaaatgttgagccttttggcctggcctgtttgcaaaggccctggccaacgtcaacattgggagcctcatctgcaa
tgtaggggcccgtggacctgctccagcagctggtgctgcaccagcaggaggtcctgccccctccactgctgctcctccagctgaggagaagaagtggaagcaagaagaagaatccgaggagctgatgatgacatggcctttggtcttttgactaa