

FKBP2-FK506 binding protein 2, 13kDa Gene

Reference: PTXBC003384

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PPIase FKBP2; peptidyl-prolyl cis-trans isomerase FKBP2; FKBP-13; PPIase; 13 kDa FK506-binding protein; 13 kDa FKBP; FK506 binding protein 2, 13kDa; FKBP-2; immunophilin FKBP13; proline isomerase; rapamycin-binding protein; rotamase; FK506 binding protein		

For research use only.

GENE DESCRIPTION	FK506 binding protein 2, 13kDa
GENE ID	2286
ORIGIN SPECIES	Human
ACCESSIONS	BC003384
VECTOR	pDONR223
NCBI SYMBOL	FKBP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggctgagctggtccgggtcctgacagtagtctccatctgcctgagcgccgtggccacggccacggggccgaggcgaaggaagctgcagatcgggtcgaagaagcgggtggaccactgtccatcaaatcgcgcaaaaggggatgctctgcacatgcactacacggggaa
gctggaagatgggacagagtttgacagcagcctgccccagaaccagcccttctctccctggcacaggccaggtcaatcaagggtgggaccaggggctgctgggatgtgtgagggggaaaagcgcaagctggtgatccatccgagctagggatggagagcggggagctcc
cccaagattccaggcgggtcaaccctggtgttcgaggtggagctgctcaaatagagcgacgaactgagctgtaa