

UBE2D3-ubiquitin-conjugating enzyme E2D 3 (UBC4/5 homolog, yeast) Gene

Reference: PTXBC003395

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E2(17)KB3; UBC4/5; UBCH5C; ubiquitin-conjugating enzyme E2 D3; (E3-independent) E2 ubiquitin-conjugating enzyme D3; E2 ubiquitin-conjugating enzyme D3; E2(17)KB 3; PRO2116; ubiquitin carrier protein D3; ubiquitin conjugating enzyme E2D 3; ubiquitin-conjug		

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2D 3 (UBC4/5 homolog, yeast)
GENE ID	7323
ORIGIN SPECIES	Human
ACCESSIONS	BC003395
VECTOR	pDONR223
NCBI SYMBOL	UBE2D3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtgaaacggattaagaacttagtgattggccgtgacccctccagcacaatgtctgcaggtccagttgggatgatattgttcattggcaagccacaattatgggacctaagacagcccatatcaaggcgggtgtattcttttgacaattcatttctacagactacccttcaa
accacctaaggttgcaattacaacaagaattatcatccaaatattaacagtaafggcagcatttgcgtatattctaagatcacagtggcgcctgtttaacaatttctaaagttctttatccatttgttactgctatgtgatccaaaccagatgacccctagtgccagagattgcaggatcta
taaacagacagagataagtacaacagaatatctcgggaatggactcagaagtatgccatgtga