



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PPIH-peptidylprolyl isomerase H (cyclophilin H) Gene

Reference: PTXBC003412

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CYP-20; CYPH; SnuCyp-20; USA-CYP; peptidyl-prolyl cis-trans isomerase H; PPIase h; U-snRNP-associated cyclophilin SnuCyp-20; U-snRNP-associated cyclophilin SunCyp-20; USA-CyP SnuCyp-20; cyclophilin H; rotamase H; small nuclear ribonucleoprotein particle-s		

For research use only.

GENE DESCRIPTION	peptidylprolyl isomerase H (cyclophilin H)
GENE ID	10465
ORIGIN SPECIES	Human
ACCESSIONS	BC003412
VECTOR	pDONR223
NCBI SYMBOL	PPIH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggtggcgaattcaagctctgttaaaccccggtgtgtcttctgattcagtaftggcggtcaggaaagtggcgcatgaagatcgagctcttgcagacgttggctcttaagacggccgagaactttaggcgattctgcaccgggaatcaggaagaatgggttccaataggatacaaa
ggaagcaccctccacaggctcataaggattcatgattcgaatggcgaggattgttttaatgagatgcttaccgagtcgcgcatttaccctggccggccatttggccagatgaaatttaaaccttcacacatcgctcagcgcctcttccatgcggcaacgttgcctcaagtaacaatggcgctgtca
ctctttcatcctctctaagctgagctgtcttgcggaaacatgtgtgttttggaaataatcatcgaatgctcttgaatgaaagaaatgagaattgtccacagagccccaacatgccaacagccactctctgtatctcagctgtggcggaatgtag

More info: Please contact us for any question: contact@proteogenix.fr