

PTMA-prothymosin, alpha Gene

Reference: PTXBC003510

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMSA; gene sequence 28; prothymosin alpha protein; prothymosin, alpha		
For research use only.			

GENE DESCRIPTION	prothymosin, alpha
GENE ID	5757
ORIGIN SPECIES	Human
ACCESSIONS	BC003510
VECTOR	pDONR223
NCBI SYMBOL	PTMA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagacgcagccgtagacaccagctccgaaatcaccaccaaggacttaaggagaagaaggaggtgtggaagagcagaaaatggaagagacgccctgctaaccgggaatgctaataagggaatggggagcaggaggctgacaatgaggtagacgaagaagaggaaga
ggggggaggaagaggaggaggaagaagaaggatggtgaggaagaggaatggagatgaagatgaggaagctgagtcagctacgggcaagcgggcagctgaagatgatgaggatgacgatgtcgataccaagaagcagaagaccgacgaggatgactag